

## Immunocytochemistry with Rabbit Polyclonal Antibodies

These procedures are for free-floating sections using trays that have bins holding approx. 5 ml. Do not fill the bins with too many sections or else staining will be uneven. Transfer with paintbrushes. All solutions are described in the Stock Solutions pdf.

### Day 1

Wash with 0.1M Tris buffer, 3 x 5min

Incubate in 1% H<sub>2</sub>O<sub>2</sub> for 1-2 min

1 ml 30% H<sub>2</sub>O<sub>2</sub> and 29 ml 0.1M Tris buffer. Always make fresh.

If bubbles form, remove the sections and put them immediately into 0.1M Tris buffer.

Wash in 0.1M Tris buffer, 3x5 min

Wash in Tris A, 10 min

Wash in Tris B, 10 min

Incubate in 10% Normal Goat Serum (S-1000, Vector Laboratories) for 30 min

4 ml Normal Goat Serum and 36 ml Tris B

Wash in Tris A, 10 min

Wash in Tris B, 10 min

Add primary antibody diluted in Tris B

Incubate overnight

### Day 2

Wash in Tris A, 10 min

Wash in Tris B, 10 min

Incubate with Anti-Rabbit IgG in Goat (secondary antibody; Vector) 1:400 for 45 min

100 µl IgG in 40 ml Tris B

Wash in Tris A, 10 min

Incubate in Tris ABC 1:1000 for 2 hrs

40 µl Reagent A

40 µl Reagent B

40 µl Tris B

Wash in Tris 3x5 min

Incubate in diaminobenzidine (DAB; Sigma)

1 DAB tablet (10 mg)

20 ml Tris buffer

40 ul Ammonium chloride stock (2g/10ml)

20 ul Glucose oxidase stock (0.0075g/2.5 ml)

160 ul D+-glucose stock (2.5g/10ml)

Wash in 0.1M Tris buffer, 3x5 min

Mount on treated slides, allow to dry overnight

Dehydrate in ethanol

70% 10 min; 90% 10 min; 95% 10 min; 100% 10 min; 100% 10 min

Clear in Xylene, 10 min

Coverslip in Permount (Fisher)

For Nickel intensification, add NiCl<sub>2</sub> to the DAB solution as follows

5 mg per 20 ml

For NovaRed (Vector) stain instead of DAB:

5 ml dH<sub>2</sub>O

3 drops Reagent 1

2 drops Reagent 2

2 drops Reagent 3

2 drops H<sub>2</sub>O<sub>2</sub> solution

wash in dH<sub>2</sub>O afterwards, 2 x 5 min, followed by 0.1M Tris buffer, 2 x 5 min