

Immunocytochemistry with Mouse Monoclonal Antibodies

This procedure is for free floating sections in trays that have bins which hold a few ml each. Make sure sections are careful transferred with paintbrushes. Limit the # sections per bin so they do not overlap and stay immersed.

Day 1

Wash with 0.1M Tris buffer (see Stock Solutions) 3 x 5 min

Incubate in 1% H₂O₂ for 1-2 min

1 ml 30% H₂O₂ and 29 ml 0.1M Tris buffer. Always make fresh.

If bubbles form, remove the sections and put them immediately into 0.1M Tris buffer

Wash in 0.1M Tris buffer 3 x 5 min

Wash in Tris A, 10 min

Wash in Tris B, 10 min

Incubate in 10% Normal Horse Serum (Vector Laboratories) for 30 min

4 ml Normal Horse Serum and 36 ml Tris B

Wash in Tris A, 10 min

Wash in Tris B, 10 min

Add primary antibody diluted in Tris B

Incubate overnight at room temperature on a rotating shaker

Day 2

Wash in Tris A, 10 min

Wash in Tris B, 10 min

Incubate with Anti-mouse IgG in Horse (secondary antibody; Vector) 1:400 for 45 min

100 µl IgG in 40 ml Tris B

Wash in Tris A, 10 min

Incubate in Tris ABC 1:1000 for 2 hrs

40 µl Reagent A

40 µl Reagent B

40 µl Tris B

Wash in 0.1M Tris buffer 3 x 5 min

Incubate in diaminobenzidine (DAB; Sigma)

1 DAB tablet (10 mg)

20 ml 0.1M Tris buffer

40 ul Ammonium chloride stock (2g/10ml)

20 ul Glucose oxidase stock (0.0075g/2.5 ml)

160 ul D+-glucose stock (2.5g/10ml)

Wash in 0.1M Tris buffer, 3 x 5 min

Mount on treated slides, allow to dry overnight

Dehydrate in ethanol sequentially:

70% 10 min; 90% 10 min; 95% 10 min; 100% 10 min; 100% 10 min

Clear in Xylene, 10 min

Coverslip in Permount (Fisher)

For Nickel intensification, add NiCl₂ to the DAB solution as follows

5 mg per 20 ml

For NovaRed (Vector) instead of DAB:

5 ml dH₂O

3 drops Reagent 1

2 drops Reagent 2

2 drops Reagent 3

2 drops H₂O₂ solution

wash in dH₂O afterwards, 2 x 5 min, followed by 0.1M Tris buffer, 2 x 5 min